

+55 (19) 2105-6161 comex@tecnal.com.br



## **OPEN CELL SHAKER** TE-083

Used to analyze disintegrated sugar cane cells.



João Leonardo Fustaino, nº 325 Distrito Industrial Uninorte Piracicaba/SP-Brasil • CEP 13.413-102



## **Technical Characteristics**

## **TE-083**

- Rotation: Fixed at 60 RPM;
- Accuracy: ± 1 RPM;
- Agitation range: 360°;
- Motor: Reducer 84 Watts;
- Capacity: 02 containers of 4 liters;
- Fixing of the containers: Through handles;
- Cabinet: In carbon steel with anti-corrosive treatment and electrostatic painting;

- Dimensions: W = 520 x D = 500 x H = 700 mm;
- Weight: 20 kg;
- Power: 100 Watts;
- Voltage: 220 Volts;
- Include: 02 Plastic containers (4 liters) 02 extra fuses Instruction Manual with Warranty Term;
- Optional: Protection cabinet with emergency stop sensor;

Benefits and Advantages

