



# ***PIPETTOR SEMI AUTOMATIC***

## ***TECSOLO-015***

Used for analysis in soil preparation for fertility purposes.

## Technical Characteristics

### TECSOLO-015

- Capillary pipettes: 11 borosilicate pipettes;
- Volume: 15 m;
- Structure: fully stainless steel;
- Dimensions: L = 918 X D=320 x H = 910 mm;
- Weight: 16kg;
- ACCOMPANY: - 45 Clamps - Warranty term;
- NOTE: it is indispensable to use vacuum pump;

### Benefits and Advantages