



SEMI-AUTOMATIC PIPETTE

TECSOLO-075

Used for analysis in soil preparation, for fertility purposes.

Technical Characteristics

TECSOLO-075

- Capillary pipettes: 11 borosilicate pipettes;
- Volume: 75 mL;
- Structure: fully stainless steel;
- Dimensions: L = 918 X D=320 x H=1015 mm;
- Weight: 16kg;
- ACCOMPANY: - 45 Clamps - Warranty term;
- Note: it is indispensable to use vacuum pump ;

Benefits and Advantages