



ROTARY SHAKER

TE-165

Used in the preparation of samples that require a 360° rotational movement.

Technical Characteristics

TE-165

- Rotation: 0 to 45 RPM;
- Accuracy: ± 3 rpm;
- Rotation control: analog;
- Shaking Amplitude: 360°;
- Motor: Reducer 38 Watts;
- Capacity: customizable proofs;
- Cabinet: in carbon steel with anti-corrosion treatment and electrostatic painting;
- Dimensions: L = 270 X D=250 x H=430 mm ;
- Weight: 7kg;
- Power: 200 Watts;
- Voltage: 220 Volts;
- ACCOMPANY: - 02 Extra fuses - Instruction Manual with Warranty Term;

Benefits and Advantages

- Promotes tube agitation without forming bubbles
- Easy operation adjusting rotation only
- Fully adjustable rotation angle
- Flexibility in choosing tubes
- Light and compact equipment
- Carbon steel cabinet with anti-corrosive treatment and electrostatic painting ensuring longer equipment life and ease of cleaning

Related Products



MAGNETIC SHAKER WITH
HEATING
TE-0854



INCUBATOR FOR BOD
TE-371/240L