



## ***EXTRACTION BATTERY***

***TE-1274***

Used for sensitivity testing of bacteria, fungi, viruses and protozoa to antimicrobials, through varying acidity.

## Technical Characteristics

### TE-1274

- Temperature: Up To 320°C;
- Temperature controller: single analog;
- Capacity: 4 proofs;
- Led: resistance Drive indicator;
- Plate: 120x120 mm aluminum;
- Resistance: in pyroceramic plate;
- Condenser Bucket: aluminum with stainless steel claws;
- Enclosure: carbon steel with anti-corrosion treatment and electrostatic painting;
- Dimensions: L = 540 x d=275 x H=775 mm;
- Power: 1600 Watts;
- Weight: 10kg;
- Voltage: 220 Volts;
- ACCOMPANY: - 04 Flat bottom balloons of 250 ml - 04 Condensers of 6 balls - Instruction Manual with Warranty Term;

## Benefits and Advantages

- Works with up to 4 Tests with independent temperature programming
- Possibility of individual use for a single sample
- Interconnected pipe system for thermostatic bath use, which ensures great water savings
- Hot areas protection system: user safety
- Quick change glassware by frosted fittings for practicality
- Good heat conservation showing a better thermal conservation from one sample to another by the porcelain resistance system
- Glass adjustment system by front and rear articulation claws for agility
- Easy maintenance
- Fast heat transfer to the sample, making the process agile
- Heating indication for analyst control and safety
- Strict quality control, in which checks and tests ensure the perfect operation of the equipment, providing safety and customer satisfaction
- Customer service, to answer questions and provide explanations about the equipment and methodologies
- Possibility of adaptations according to the needs of the customer, which makes the equipment already in line a special equipment.