



SOIL PHYSICAL ANALYSIS STIRRER

TE-161

Rotary mechanical stirrer used to stir samples for subsequent physical soil/granulometric analysis.

Technical Characteristics

TE-161

- Rotation: 0 to 60 RPM;
- Accuracy: ± 5 RPM;
- Rotation control: Analog with digital display;
- Agitation range: 360°;
- Engine: Induction 1/3 HP;
- Timer: Digital - Programmable up to 99:59 hours. Automatic shutdown at the end of the scheduled time;
- Capacity: 60 tests;
- Support: Swivel with 10 (ten) drawers for 500 ml test tubes;
- Fixing the drawers: Through handles;
- Cabinet: In carbon steel with anti-corrosive treatment and electrostatic painting;
- Dimensions: W=1150 x D=650 x H=670 mm;
- Weight: 80kg (without sample);
- Power: 250 Watts;
- Voltage: 220 Volts;
- ACCOMPANIES: - 02 extra fuses - 60 500ml beakers - Instructions Manual and Warranty Term;

Benefits and Advantages

- Acrylic protection, ensuring user safety
- Supports large number of samples, optimizing the analysis process
- Emergency button for safety, which stops the operation of the equipment when activated
- Automatic shutdown at the end of the programmed time
- Resistant material, ensuring a longer useful life for the equipment
- Rigid Quality Control, in which checks and tests guarantee the perfect functioning of the equipment, providing safety and client satisfaction.

Related Products



WAGNER SHAKER
TE-160/24



VERTICAL TEST TUBE
STIRRER
TE-167